

# DEPT Experiment Guide

Step	Function or Dialog Box	<Keystroke>/[Select]<Data Entry>	Comment
1	Sample		See Sample Preparation Guide. Position sample spinner using the depth gauge, place in probe.
2	Enter PNMR program.	<Alt+Tab>	(If necessary.)
3	Select $^{13}\text{C}$ observe.	<b>H1&gt;nu C13&lt;Enter&gt;</b>	Required only if the prompt is not <b>C13&gt;</b> .
	<i>Optional:</i> Shim sample.	<b>C13&gt;shim&lt;Enter&gt;</b>	Follow on-screen directions. Shims are optimized.
4	Acquire $^1\text{H}$ spectrum.	<b>C13&gt;zgh&lt;Enter&gt;</b>	
5	Enter NUTS and process data.	<Alt+Tab> >a2	Process $^1\text{H}$ spectrum with a2 link. Trim phase as required. Use cursor to determine TMS peak position in ppm, including sign.
6	Enter PNMR program.	<Alt+Tab>	
7	Enter TMS peak position.	<b>C13&gt;fo&lt;Enter&gt;</b> <b>value&lt;Enter&gt;</b> <b>0&lt;Enter&gt;</b>	In the first dialog box enter the current position (in ppm) of the TMS peak. In the second dialog box enter 0 (zero). <b>Repeat to confirm.</b>
8	Verify parameters.		Verify that parameters make sense; for neat samples NS=12, for 1M samples NS=60.
9	Run DEPT program.	<b>C13&gt;dept&lt;Enter&gt;</b> then <i>filename&lt;Enter&gt;</i> or <Enter> for default	The program requests a file name while providing the default name My_dept. Use the default unless planning to save the data long term.
10	Acquire data.		The pulse program runs DEPT45, DEPT90, and DEPT135, requiring about 3 min using the parameters specified in step 8.
11	Enter NUTS.	<Alt+Tab>	
12	Process data.	<Ctrl+F11> then [filename][Open] to select a file or [Open] for default	Data is processed with macro aii_dept.mac. Use mouse to select default data file [my_dept] or named data file [filename].
13	Enter line broadening.	<b>value&lt;Enter&gt;</b>	<b>LB = 0.5 Hz is a typical value.</b>
14	Display Stacked Plot.	↑ and/or ↓	Macro applies an approximate chemical shift scale, sets display for 220 to -10 ppm window, and displays a stacked plot of the DEPT45, DEPT90, and DEPT135 spectra. Adjust amplitude with Page Up/Page Down, arrow keys or vertical slider
15	Print Stacked Plot.	<p>	
16	Exit Stacked Plot.	<Enter>	To redisplay stacked plot from the base level NUTS prompt, use the command “sp”.