

Table of Contents

3	2a	Measuring small volumes in a Biotechnology Lab
5	2b	measuring very small volumes in a Biotechnology Lab
9	2c	Measuring Mass
11	2d	Checking the Accuracy of the Micropipets using a balance
13	3a	Making solutions of differing mass/volume concentrations
15	3b	Making solutions of differing % mass/volume concentrations
17	3c	Making solutions of differing molarity concentrations
19	3d	Making dilutions of concentrated solutions
23	4a	Measuring the pH of solutions.
25	4b	Making a buffer solution at a specific pH
27	4c	Demonstration of Buffer Efficacy
33	4d	DNA isolation from Strawberries.
35	5a	Restriction Enzyme Digestion of DNA
37	5b	Gel Electrophoresis
39	5c	DNA staining
45	6a	Assaying for Starch and Sugar
49	6b	Assaying for Amylase activity
51	6c	Assaying Plant Substances for Antibacterial Activity
53	7a	Data Analysis and Conclusions: Amylase Activity Assay
67	7b	Assaying antibacterial activity in plant extracts
61	7c	Assaying different enzymes for apple juice production
63	8a	Testing samples for Hydrogen Peroxidase activity
65	8b	Isolation of HRP from horseradish root.
67	8c	Testing for the presence of Peroxidase using TMB
69	9a	Bacterial transformation
71	9b	using spectrophotometer
73	9c	using spectrophotometer to study molecules
77	10a	Transformation, calculating transformation efficiency
81	11a	GFP purification
83	11b	Using the spectrophotometer to study protein concentration
87	11c	Determining the concentration of BSA Protein in solution
91	12a	GFP Purification
93	12b	Determination of Protein concentration
97	12c	Preparation of samples for SDS PAGE
101	13a	SDS-PAGE sample preparation and Gel loading
103	13b	Protein gel staining & western blot
105	14a	Western Blot Part 2: color development
109	14b	Analysis of results from protein gel staining

Title: 2a Measuring small volumes in a Biotechnology Lab

Reference: 14-15 lab manual

Purpose: To learn how to measure small volumes correctly with pipets.

Procedure: Label each of your plastic tubes A, B, C + D. Make sure to understand the maximum amount of liquid each pipet can measure. Measure the correct amount of each dye into the tubes A, B, C, D and check that your tubes are within the meniscus of the key tube. Redo if not.

Tube	Red dye vol (mL)	Blue dye vol (mL)	Yellow dye vol (mL)	Total volume (mL)
A	1.1	1.7	0.33	3.1
B	1.27	0.85	2.9	5.0
C	0.7	2.8	1.8	5.3 5.0
D	1.6	1.9	0.66	

Observations and Procedural Notes

It took some practice to figure out which pipette to use. The right size for accuracy and ease of use. We used the 1-2 mL for this lab with .01 mL graduations. And for one the 5 mL for one. After practice we all got much more efficient w/ using the electric pipettor. Speed of suction was crucial. After a few tries we all learned to use a gentle, practiced hand. As time went on we were able to pipette much slower. One was close but I added a mL by starting at -1

Title: Measuring very small volumes in a biotechnology Lab 26

Reference: Pages 15-17 Lab manual

Procedure:

Purpose: to learn how to measure very small volumes in μL with micropipettes and then calculate air percent error.

Procedure: Following Table 2.2 and using the correct micropipette for each amount, measure the indicated amount of dye into each tube. Spin the tubes in the microcentrifuge. Measure the amount of liquid you pipetted into each tube and calculate your percent error.

Table 2.3 Percent error for micropipetting practice

Tube	Expected volume	Actual volume	% Error
A	12 μL	13 μL	8.3%
B	24 μL	21 μL	12.5%
C	100 μL	105	5.0%
D	900 μL	855	5.0%

13%

Observations and Procedural notes

The micropipette was much easier to use once we mastered the first & second stop.

With Tube A with measuring the actual volume we had some liquid left with the micropipette set at 12 μL so we switched to 13 μL and then had air bubbles so we made an error.

Tube error was measured

6

ONCE we ~~measured~~
mastered measuring the
volume from the micropipette
tubes we were able to
measure the actual amount we
measured much faster.

Hopefully we can improve
our accuracy with time.

$$\text{tube A} \quad \left| \frac{a - e}{e} \right| \times 100 \text{ error}$$

$$\left| \frac{13 - 12}{12} \right| \times 100 = 8.3\%$$

$$\text{tube C} \quad \left| \frac{105 - 100}{100} \right| \times 100 = 5.0\%$$

$$\text{tube D} \quad \left| \frac{955 - 900}{900} \right| \times 100 = 5.0\%$$

$$\text{tube B} \quad \left| \frac{24 - 24}{24} \right| \times 100 = 12.5\% \\ = 13\%$$

Title: Measuring Mass 2c

Reference: Lab manual pages 17-19

Purpose: To learn how to accurately measure mass using an analytical scale

Procedure: weigh out the appropriate amount of glucose according to the chart for each solution. Add water to each tube to make a final volume of 10 mL. Test each tube with a diastix test strip and compare color to the key. Record the results.

Solution to be prepared mg/dL	Equivalent concentration g/mL	Final volume (mL)	mass of glucose (g)	Diastix test result
250.0	0.0025	10	0.025	225 mg/dL
100.0	0.001	10	0.01	100 mg/dL
50.0	0.0005	10	0.005	50 mg/dL
0	0	10	0	0 mg/dL

Observations and Procedural notes

We had a little trouble with the analytical balance at first, until we realized it took a very light touch. Weigh paper would normally work better but we measured with the small weigh boats. The only sample where we were not sure was the 250.0 mg/dL it was slightly off so we gave it a value of 225 mg/dL. It could have

Title: 2.2 Checking the Accuracy of Micropipets using a Balance

Reference: Lab Manual 18-19

Purpose: To measure the accuracy of micropipets by using the known factor of 1ml water weighs 1.0g

Procedure: measure the specified volume of water into the weigh boats and then weigh the water. calculate the percent error. remeasure if outside of accepted range.

Table 2.5 micropipetting precision

Micropipet	Volume (µL)	Volume (mL)	Expected Mass (g)	Actual mass (g)	Percent Error	Acceptable error (%)	Type of Balance used
P-1000	1000	1.0	1.0	1.03	3%	3	table top
P-1000	255	.255	.255	.278 .2508	1.65%	3	analytical
P-100	100	1.00	.100	.100	0%	5	analytical
P-100	17	.017	.017	.018 .0178	4.71%	5	analytical
P-10	10	.01	.01	.011	10%	10	analytical
P-10	2.4	.0024	.0024	.003 .0025	4%	10	analytical

$$\left| \frac{\text{actual mass} - \text{expected mass}}{\text{expected mass}} \right| \times 100$$

$$\left| \frac{1.03 - 1.0}{1.0} \right| \times 100 = 3\%$$

the first P-1000 required recalibration after we were outside the recommended error $> 10\%$

The next P-1000 was accurate.

- the P-100 was more accurate than the P-1000 for measuring 100 μL

if the measurement ~~if~~ was outside of acceptable error, we redid the ~~calculation~~ measurement.

$$\left| \frac{.100 - .100}{.100} \right| \times 100 = 0\%$$

$$\left| \frac{.0178 - .017}{.017} \right| \times 100 = 4.7\% \quad 0\%$$

$$\left| \frac{.011 - .01}{.01} \right| \times 100 = 10\%$$